Roche analyser / module guidance for registration in UK NEQAS Clinical Chemistry Schemes

Doc Code: S-WKS-1138
Version: 2

This guidance has been written to assist with Roche method classification for the Birmingham Quality, UK NEQAS Cardiac Markers and UK NEQAS Guildford Peptides partnership. It has ***not*** been endorsed by Roche. It should ***only*** be used for External Quality Assessment (EQA) purposes.

We recognise that there are multiple Roche systems available:

- cobas® 6000 series — c501 *and/or* e601
- cobas® 8000 series — c502/c701/c702 *and/or* e602/e801
- cobas® pure integrated solutions — c303 *and/or* e402
- cobas® pro integrated solutions — c503/c703 *and/or* e801
- c111 — clinical chemistry only
- c311 — clinical chemistry only
- e411 — immunoassay only

For EQA purposes, we will apply the following method classifications in the new partnership Schemes.

Clinical Chemistry Analyser/Module

- Roche URC I (311/501/502)
- Roche URC II (701/702)
- Roche URC III (303/503/703)
- Roche reagent bottles (111)

Immunoassay Analyser/Module

- Roche E1G (411/601/602)
- Roche E2G (402/801)

If you have any questions, please contact birminghamquality@uhb.nhs.uk